


Transpulmonary thermodilution TPTD



For professional use. CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

Edwards Lifesciences devices placed on the European market meeting the essential requirements referred to in Article 3 of the Medical Device Directive93/42/ECC bear the CE marking of conformity.

Edwards, Edwards Lifesciences, and the stylized E logo are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

© 2017 Edwards Lifesciences Corporation.
All rights reserved. PP--US-2312 v1.0

Edwards Lifesciences
One Edwards Way, Irvine CA 92614 USA
edwards.com



Normal Hemodynamic Parameters and Laboratory Values

Normal hemodynamic parameters – adult

Parameter	Equation	Normal range
Arterial Oxygen Saturation (SaO ₂)		95-100%
Mixed Venous Saturation (SvO ₂)		60-80%
Central Venous Oxygen Saturation (ScvO ₂)		60-70%
Arterial Blood Pressure (BP)	Systolic (SBP) Diastolic (DBP)	90-140 mmHg 60-90 mmHg
Mean Arterial Pressure (MAP)	SBP + (2 x DBP)/3	70-105 mmHg
Right Atrial Pressure (RAP)		2-6 mmHg
Central Venous Pressure (CVP)		2-6 mmHg
Right Ventricular Pressure (RVP)	Systolic (RVSP) Diastolic (RVDP)	15-25 mmHg 0-8 mmHg
Pulmonary Artery Pressure (PAP)	Systolic (PASP) Diastolic (PADP)	15-25 mmHg 8-15 mmHg
Mean Pulmonary Artery Pressure (MPAP)	PASP + (2 x PADP)/3	10-20 mmHg
Pulmonary Artery Occlusion Pressure (PAOP)		6-12 mmHg
Left Atrial Pressure (LAP)		6-12 mmHg
Cardiac Output (CO)	HR x SV/1000	4.0-8.0 L/min
Cardiac Index (CI)	CO/BSA	2.5-4.0 L/min/m ²
Stroke Volume (SV)	CO/HR x 1000	60-100 mL/beat
Stroke Volume Index (SVI)	CI/HR x 1000	33-47 mL/m ² /beat
Stroke Volume Variation (SVV)	(SVmax - SVmin)/SVmean x 100	10-15%
Systemic Vascular Resistance (SVR)	80 x (MAP - RAP)/CO	800-1200 dynes-sec/cm ⁵
Systemic Vascular Resistance Index (SVRI)	80 x (MAP - RAP)/CI	1970-2390 dynes-sec/cm ⁵ /m ²
Pulmonary Vascular Resistance (PVR)	80 x (MPAP - PAOP)/CO	<250 dynes-sec/cm ⁵
Pulmonary Vascular Resistance Index (PVRI)	80 x (MPAP - PAOP)/CI	255-285 dynes-sec/cm ⁵ /m ²
Left Ventricular Stroke Work Index (LVSWI)	SVI x (MAP - PAOP) x 0.0136	50-62 mmHg x mL/m ²
Right Ventricular Stroke Work Index (RVSWI)	SVI x (MPAP - CVP) x 0.0136	5-10 mmHg x mL/m ²



Normal hemodynamic parameters – adult

Parameter	Equation	Normal range
Coronary Artery Perfusion Pressure (CPP)	Diastolic BP-PAOP	60-80 mmHg
Right Ventricular End-Diastolic Volume (RVEDV)	SV/EF	100-160 mL
Right Ventricular End-Diastolic Volume Index (RVEDVI)	RVEDV/BSA	59-94 mL/m ²
Right Ventricular End-Systolic Volume (RVESV)	EDV-SV	50-100 mL
Right Ventricular Ejection Fraction (RVEF)	SV/EDV x 100	40-60%
Arterial Oxygen Content (CaO ₂)	(0.0138 x Hgb x SaO ₂) + 0.0031 x PaO ₂	17-20 mL/dL
Venous Oxygen Content (CvO ₂)	(0.0138 x Hgb x SvO ₂) + 0.0031 x PvO ₂	12-15 mL/dL
A-V Oxygen Content Difference (C(a-v)O ₂)	CaO ₂ - CvO ₂	4-6 mL/dL
Oxygen Delivery (DO ₂)	CaO ₂ x CO x 10	950-1150 mL/min
Oxygen Delivery Index (DO ₂ I)	CaO ₂ x CI x 10	500-600 mL/min/m ²
Oxygen Consumption (VO ₂)	C(a-v)O ₂ x CO x 10	200-250 mL/min
Oxygen Consumption Index (VO ₂ I)	C(a-v)O ₂ x CI x 10	120-160 mL/min/m ²
Oxygen Extraction Ratio (O ₂ ER)	(CaO ₂ - CvO ₂)/CaO ₂ x 100	22-30%
Oxygen Extraction Index (O ₂ EI)	(SaO ₂ - SvO ₂)/SaO ₂ x 100	20-25%

Normal blood laboratory values

Test	Convention units (reference values*)	SI units
Hematocrit (Hct)	Males: 42-52% Females: 36-48%	0.42-0.52 0.36-0.48
Hemoglobin (Hgb)	Males: 12.4-17.4 g/dL Females: 11.7-16 g/dL	124-174 g/L 117-160 g/L
Lactate	0.93-1.65 mEq/L	0.93-1.65 mmol/L

SI Units = International Units

*Reference Values vary by regional laboratory techniques and methods.

Normal hemodynamic parameters – adult

Parameter	Equation	Normal range
Extra Vascular Lung Water (EVLW)	CO x DSt - 0.25xGEDV	
Extra Vascular Lung Water Index (ELWI)	EVLW/PBW Predicted Body Weight (PBW): Female: 45.5 + 0.91 x (Height-152.4) Male: 50 + 0.91 x (Height-152.4)	0-7 mL/kg
Global End Diastolic Volume (GEDV)	CO x MTt x f(S1/S2)	
Global End Diastolic Volume Index (GEDI)	CO x MTt x f(S1/S2)	650-800 mL/kg
Global Ejection Fraction (GEF)	SV x 4 / GEDV	>20%
Cardiac Function Index (CFI)	1000 x CO / GEDV	4.5-6.5 1/min
Intra Thoracic Blood Volume (ITBV)	ITBV = 1.25 x GEDV	
Intra Thoracic Blood Volume Index (ITBI)	ITBI = 1.25 x GEDI	850-1000 mL/m ²
Pulmonary Vascular Permeability Index (PVPI)	EVLW/0.25 x GEDV	<3
Cardiac Power (CPO)	CO x MAP x K	
Cardiac Power Index (CPI)	CI x MAP x K	0.5-0.7 W/m ²

Physio-relationship

