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
- Akkermans J, Diepeveen M, Ganzevoort W, van Montfrans GA, Westerhof BE, Wolf H. Continuous non-invasive blood pressure monitoring, a validation study of Nexfin in a pregnant population. *Hypertens Pregnancy*. 2009 May; 28(2):230-42. doi: 10.1080/10641950802601260.
- Eeftinck Schattenkerk DW, Van Lieshout JJ, Van den Meiracker AH, Wesseling KR, Blanc S, Wieling W, Van Montfrans GA, Settels JJ, Wesseling KH, Westerhof BE. Nexfin noninvasive continuous blood pressure validated against Riva-Rocci/Korotkoff. *Am J Hypertens*. 2009;22:378-83.
- Martina JR, Westerhof BE, Van Goudoever J, de Beaumont EM, Truijen J, Kim YS, Immink RV, Jobsis DA, Hollmann MW, Lahpor JR, De Mol BA, Van Lieshout JJ. Noninvasive continuous arterial blood pressure monitoring with nexfin. *Anesthesiology*. 2012;116:1092-103.
- Fischer MO, Avram R, Cărijaliu I, Massetti M, Gérard JL, Hanouz JL, Fellahi JL. Non-invasive continuous arterial pressure and cardiac index monitoring with Nexfin after cardiac surgery. *Br J Anaesth*. 2012 Oct;109(4):514-21.
- Martina JR, Westerhof BE, Van Goudoever J, De JN, Van Lieshout JJ, Lahpor JR, De Mol BA. Noninvasive blood pressure measurement by the Nexfin monitor during reduced arterial pulsatility: a feasibility study. *ASAIO J*. 2010;56:221-7.
- Vos JJ, Poterman M, Mooyaart EA, Weening M, Struys MM, Scheeren TW, Kalmars AF. Comparison of continuous non-invasive finger arterial pressure monitoring with conventional intermittent automated arm arterial pressure measurement in patients under general anaesthesia. *Br J Anaesth*. 2014 Jul;113(1):67-74.
- Bubeneck-Turconi SI, Craciun M, Miclea I, Perel A. Noninvasive Continuous Cardiac Output by the Nexfin Before and After Preload-Modifying Maneuvers: A Comparison with Intermittent Thermodilution Cardiac Output. *Anesth Analg*. 2013 Aug;117(2):366-72.
- Bogert LW, Wesseling KH, Schraa O, Van Lieshout EJ, De Mol BA, Van GJ, Westerhof BE, Van Lieshout JJ. Pulse contour cardiac output derived from non-invasive arterial pressure in cardiovascular disease. *Anaesthesia*. 2010;65:1119-25.
- Stover JF, Stocker R, Lenherr R, Neff TA, Cottini SR, Zoller B, Béchir M. Noninvasive cardiac output and blood pressure monitoring cannot replace an invasive monitoring system in critically ill patients. *BMC Anesthesiol*. 2009 Oct 12;9:6. doi: 10.1186/1471-2253-9-6.
- Broch O, Renner J, Gruenewald M, Meybohm P, Schöttler J, Caliebe A, Steinfath M, Malbrain M, Bein B. A comparison of the Nexfin and transcardiopulmonary thermodilution to estimate cardiac output during coronary artery surgery. *Anaesthesia* 2012;67:377-83.
- Hofhuizen C, Lansdorp B, van der Hoeven JG, Scheffer GJ, Lemson J. Validation of noninvasive pulse contour cardiac output using finger arterial pressure in cardiac surgery patients requiring fluid therapy. *Journal of Critical Care* 2014 Feb;29(1):161-5.
- van der Spoel AG, Voogel AJ, Folkers A, Boer C, Bouwman RA. Comparison of noninvasive continuous arterial waveform analysis (Nexfin) with transthoracic Doppler echocardiography for monitoring of cardiac output. *J Clin Anesth*. 2011 Jun;24(4):304-9. doi: 10.1016/j.jclinane.2011.09.008.
- Chen G, Meng L, Alexander B, Tran NP, Kain ZN, Cannesson M. Comparison of noninvasive cardiac output measurements using the Nexfin monitoring device and the esophageal Doppler. *J Clin Anesth*. 2012 Jun;24(4):275-83. doi: 10.1016/j.jclinane.2011.08.014.
- Felix Balzer, Marit Habicher, Michael Sander, Julian Sterr, Stephanie Scholz, Aarne Feldheiser, Michael Müller, Carsten Perka and Sascha Treskatsch. Comparison of the noninvasive Nexfin monitor with conventional methods for the measurement of arterial blood pressure in moderate risk orthopedic surgery patients. *Journal of International Medical Research*, 44(4):832-43, 2016.
- R.B.P. de Wilde, F. de Wit, B.F. Geerts, A.L. van Vliet, L.P.H. Aarts, J. Vuyk and J.R.C. Jansen. Non-invasive continuous arterial pressure and pulse pressure variation measured with Nexfin in patients following major upper abdominal surgery: a comparative study. *Anaesthesia* 2016, 71, 788-797.
- J.F. Heusdens, S. Lof, C.W.A. Pennekamp, J.C. Specken-Welleweerd, G.J. de Borst, W.A. van Klei, L. van Wolfswinkel and R.V. Immink. Validation of non-invasive arterial pressure monitoring during carotid endarterectomy. *British Journal of Anaesthesia*, 117 (3): 316-23 (2016).
- G. F. N. Berkelmans, S. Kuipers, B. E. Westerhof, A. M. E. Spoelstra, de Man, Y. M. Smulders. Comparing volume-clamp method and intra-arterial blood pressure measurements in patients with atrial fibrillation admitted to the intensive or medium care unit. *J Clin Monit Comput* DOI 10.1007/s10877-017-0044-9.
- Dorothea E. Rogge, Julia Y. Nicklas, Gerhard Schön, Oliver Grothe, Sebastian A. Haas, Daniel A. Reuter, and Bernd Saugel. Continuous Noninvasive Arterial Pressure Monitoring in Obese Patients During Bariatric Surgery: An Evaluation of the Vascular Unloading Technique (ClearSight system). 2019 Mar;128(3):477-483. doi: 10.1213/ANE.0000000000003943.
- Alberto Noto, Filippo Sanfilippo, Giovanni De Salvo, Claudia Crimi, Filippo Benedetto, Ximena Watson, Maurizio Ceconi and Antonio David. Noninvasive continuous arterial pressure monitoring with ClearSight during awake carotid endarterectomy. *Eur J Anaesthesiol* 2018; 35:1-9.

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Edwards Lifesciences • Route de l'Étraz 70, 1260 Nyon, Switzerland • edwards.com

Edwards

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© 2020 Edwards Lifesciences Corporation. Todos los derechos reservados. PP--EU-0377 v1.0