




Experience the Difference



Continuous ScvO₂ monitoring is an early indicator and key to assessing the adequacy of the balance of oxygen delivery and consumption. The goal of continuous ScvO₂ monitoring with the PreSep oximetry catheter is to bring into balance the relationship between oxygen delivery and consumption to improve the care of high-acuity patients. The PreSep oximetry catheter kit has been updated with two new components.

Experience the Difference...






Nitinol Guidewire

Designed to ease insertion and reduce extraction force.


- Kink resistant
- Flexible tip



Dilator

Designed to reduce insertion force providing greater control and reducing incision size.

- Smaller inner diameter (0.032")
- Coated
- Reduced taper angle



For professional use. CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

Edwards Lifesciences devices placed on the European market meeting the essential requirements referred to in Article 3 of the Medical Device Directive 93/42/EEC bear the CE marking of conformity.

Edwards, Edwards Lifesciences, the stylized E logo, and PreSep are trademarks of Edwards Lifesciences Corporation.

© 2012 Edwards Lifesciences Corporation. All rights reserved. AR08011

Edwards Lifesciences | edwards.com

Irvine, USA | Nyon, Switzerland | Tokyo, Japan | Singapore, Singapore | São Paulo, Brazil

