

Fogarty® Occlusion Catheters

Temporary vessel occlusion

Edwards Lifesciences offers the Fogarty occlusion catheter in a variety of sizes.



Edwards

Model number	Diameter of inflated balloon	Catheter French size	Maximum French size over deflated balloon	Length
620403F	5 mm	3F	4.3F	40 cm
620404F	9 mm	4F	5.0F	40 cm
620405F	11 mm	5F	6.0F	40 cm
62080814F	28 mm	8F	14F	80 cm
62080822F	45 mm	8F	22F	80 cm

When you are seeking time-tested solutions for temporary occlusion, reach for Fogarty occlusion catheters.

CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

Edwards, Edwards Lifesciences, the stylized E logo, and Fogarty are trademarks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2020 Edwards Lifesciences Corporation. All rights reserved. PP--US-0660 v2.0

Edwards Lifesciences • One Edwards Way, Irvine CA 92614 USA • edwards.com



Edwards